



Product Datasheet

U2SURP antibody (orb390288)


Description	Rabbit polyclonal antibody to U2SURP.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human U2SURP(NP_001073884.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	140kDa
Uniprot ID	O15042
Dilution Range	WB: 1:500-2000, IHC-P: 1:20-200, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse li...



Immunohistochemical staining of human li...



Immunohistochemical staining of rat live...