



Product Datasheet

SKIL antibody (orb390300)

Description	Rabbit polyclonal antibody to SKIL.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human SKIL(NP_005405.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75KDa
Uniprot ID	P12757
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human skeletal muscle tissue.



Western blot analysis of extracts of various tissues.