



Product Datasheet

ATG7 antibody (orb390306)


Description	Rabbit polyclonal antibody to ATG7.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human ATG7(NP_001129503.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68kDa/75kDa/77kDa
Uniprot ID	O95352
Dilution Range	IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric cancer tissue showing ATG7 expression.



Immunohistochemical staining of human breast cancer tissue showing ATG7 expression.



Immunohistochemical staining of human colon cancer tissue showing ATG7 expression.