



Product Datasheet

EEF1G antibody (orb390311)


Description	Rabbit polyclonal antibody to EEF1G.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human eEF1G(NP_001395.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50kDa
Uniprot ID	P26641
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH-3T3 cells



Immunofluorescence analysis of C6 cells



Immunohistochemical staining of human esophageal tissue