


Product Datasheet

ZRANB3 antibody (orb390328)

Description	Rabbit polyclonal antibody to ZRANB3.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ZRANB3(NP_001273497.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	123kDa
Uniprot ID	Q5FWF4
Dilution Range	WB: 1:1000-3000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of various cell lines.