



Product Datasheet

RELA (Phospho-S311) antibody (orb390336)

Description	Rabbit polyclonal antibody to RELA (Phospho-S311).
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S311 of human RelA(NP_068810.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	65KDa
Uniprot ID	Q04206
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse sp...



Western blot analysis of extracts of HeL...