

Product Datasheet

Phenylalanine ammonia-lyase Microplate Assay Kit (orb390724)

Description	Detection and Quantification of Phenylalanine ammonia-lyase Activity.
Reactivity	All
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: ColorimetricSample Types: Tissue extracts, Cell lysate
Sample Types	Tissue extracts, Cell lysate
CAS Number	9024-28-6
Expiration Date	12 months from date of receipt.