



Product Datasheet

Alpha-Amylase Microplate Assay Kit (orb390729)

Description	Detection and Quantification of Alpha-Amylase Activity.
Reactivity	All
Range	0.2 mmol/L - 2 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate
Sample Types	Tissue extracts, Cell lysate
CAS Number	9000-90-2
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...