

Product Datasheet

Beta-Amylase Microplate Assay Kit (orb390730)

Description Detection and Quantification of Beta-Amylase Activity.

Reactivity All

Range 0.2 mmol/L - 2 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate


Sample Types Tissue extracts, Cell lysate

CAS Number 9000-91-3

Expiration Date 12 months from date of receipt.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...