


Product Datasheet

Glutathione S-transferase Microplate Assay Kit (orb390753)

Description	Detection and Quantification of Glutathione S-transferase Activity.
Reactivity	All
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	50812-37-8
Expiration Date	12 months from date of receipt.



Positive Control reaction in 96-well pla...