


Product Datasheet

Ascorbic Acid Microplate Assay Kit (orb390754)

Description	Detection and Quantification of Ascorbic Acid Content.
Reactivity	All
Range	0.02 mmol/L - 2 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate
Sample Types	Tissue extracts, Cell lysate
CAS Number	50-81-7
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve