


Product Datasheet

Dehydroascorbic Acid Microplate Assay Kit (orb390755)

Description	Detection and Quantification of Dehydroascorbic Acid Content.
Reactivity	All
Range	50 µmol/L - 500 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate
Sample Types	Tissue extracts, Cell lysate
CAS Number	490-83-5
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve