

Product Datasheet

L-galactono-1,4-lactone Dehydrogenase Microplate Assay Kit (orb390756)

Description	Detection and Quantification of L-galactono-1,4-lactone Dehydrogenase Activity.
Reactivity	All
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: ColorimetricSample Types: Tissue extracts, Cell lysate
Sample Types	Tissue extracts, Cell lysate
Expiration Date	12 months from date of receipt.