

Product Datasheet

PIK3CG antibody (orb39076)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	PIK3CG
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	126454
Uniprot ID	P48736
NCBI	NP_002640.2, NP_001269356.1, NP_001269355.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

