

Product Datasheet

Diamine Oxidase Microplate Assay Kit (orb390761)

Description Detection and Quantification of Dehydroascorbate Reductase Activity.

Reactivity All

Range 0.05 mmol/L - 5 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids


Sample Types Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

CAS Number 9001-53-0

Expiration Date 12 months from date of receipt.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...