

Product Datasheet

Glucose Oxidase Microplate Assay Kit (orb390763)

Description Detection and Quantification of Glucose Oxidase Activity.

Reactivity All

Range 0.03 mmol/L - 3 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids

Sample Types Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

CAS Number 9001-37-0

Expiration Date 12 months from date of receipt.



Representative
Standard
Curve