


Product Datasheet

Uric Acid Microplate Assay Kit (orb390764)

Description	Detection and Quantification of Uric Acid Content.
Reactivity	All
Range	0.04 mmol/L - 4 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Urine, Serum, Plasma, Other biological fluids
Sample Types	Urine, Serum, Plasma, Other biological fluids
CAS Number	69-93-2
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve