

Product Datasheet

Peroxidase Microplate Assay Kit (orb390768)

Description Detection and Quantification of Peroxidase Activity.

Reactivity All

Tested Applications FA


Storage Shipped and store at 4°C for 6 months.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids

Sample Types Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

Expiration Date 12 months from date of receipt.



Positive Control reaction in 96-well pla...