



Product Datasheet

Glutamate Dehydrogenase Microplate Assay Kit (orb390772)

Description	Detection and Quantification of Glutamate Dehydrogenase Activity.
Reactivity	All
Range	4 µmol/L - 400 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate
Sample Types	Tissue extracts, Cell lysate
CAS Number	9001-46-1
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...