


Product Datasheet

Hydroxyproline Microplate Assay Kit (orb390773)

Description	Detection and Quantification of Hydroxyproline Content.
Reactivity	All
Range	2 µmol/L - 200 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate
Sample Types	Tissue extracts, Cell lysate
CAS Number	51-35-4
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve