



Product Datasheet

Acetylcholinesterase Microplate Assay Kit (orb390774)

Description	Detection and Quantification of Acetylcholinesterase Activity.
Reactivity	All
Range	0.1 mmol/L - 5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9000-81-1
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...