

Product Datasheet

Citrate Microplate Assay Kit (orb390775)

Description Detection and Quantification of Citrate Content.

Reactivity All

Range 0.05 mmol/L - 5 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids

Sample Types Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

CAS Number 77-92-9

Expiration Date 12 months from date of receipt.



Positive
Control
reaction in
96-well pla...