

Product Datasheet

Alpha-Glucosidase Microplate Assay Kit (orb390802)

Description Detection and Quantification of Alpha-Glucosidase Activity.

Reactivity All

Range 10 µmol/L - 1000 µmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Tissue extracts, Cell lysate, Cell culture media and Other biological fluids


Sample Types Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

CAS Number 9001-42-7

Expiration Date 12 months from date of receipt.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...