



Product Datasheet

Beta-Glucosidase Microplate Assay Kit (orb390803)

Description	Detection and Quantification of Beta-Glucosidase Activity.
Reactivity	All
Range	10 µmol/L - 1000 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-22-3
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...