


Product Datasheet

Monoamine Oxidase Microplate Assay Kit (orb390808)

Description	Detection and Quantification of Monoamine Oxidase Activity.
Reactivity	All
Range	0.1 mmol/L - 10 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-66-5
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve