

Product Datasheet

Serum Kalium Microplate Assay Kit (orb390810)

Description Detection and Quantification of Serum Kalium Content.

Reactivity All

Range 0.2 mmol/L - 2 mmol/L

Tested Applications FA


Storage Shipped and store at 4°C for 6 months.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Serum

Sample Types Serum, Urine, Saliva, Other biological fluids

Expiration Date 12 months from date of receipt.



Representative
Standard
Curve