


Product Datasheet

Serum Phosphorus Microplate Assay Kit (orb390814)

Description	Detection and Quantification of Serum Phosphorus Content.
Reactivity	All
Range	0.01 mmol/L - 0.4 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Serum
Sample Types	Serum, Urine, Saliva, Other biological fluids
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve