


Product Datasheet

Total Antioxidant Capacity Microplate Assay Kit (orb390817)

Description	Detection and Quantification of Total Antioxidant Capacity Activity.
Reactivity	All
Range	0.04 mmol/L - 4 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve