

## Product Datasheet

### Human ADIPOQ ELISA Kit (orb390821)

**Description** Sandwich ELISA for quantitative detection of human ADIPOQ concentrations.

**Reactivity** Human

**Range** 62.5 pg/ml - 4000 pg/ml

**Tested Applications** ELISA

**Target** ADIPOQ

**Storage** Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

**Note** For research use only

**Application notes** Detection Method: Colorimetric  
Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates


**Sample Types** Cell Culture Supernates, Serum, Plasma, Tissue Homogenates

**Uniprot ID** **Q15848**

**Sensitivity** 15 pg/ml

**Entrez** **9370**

**Expiration Date** Please enquire.



Representative  
Standard  
Curve