

Product Datasheet

Human ADIPOQ ELISA Kit (orb390822)

Description Sandwich ELISA for quantitative detection of human ADIPOQ concentrations.

Reactivity Human

Range 62.5 pg/ml - 4000 pg/ml

Tested Applications ELISA

Target ADIPOQ

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates


Sample Types Cell Culture Supernates, Serum, Plasma, Tissue Homogenates

Uniprot ID **Q15848**

Sensitivity 15 pg/ml

Entrez **9370**

Expiration Date Please enquire.



Representative
Standard
Curve