

Product Datasheet

Human ANGPT1 ELISA Kit (orb390824)

Description Sandwich ELISA for quantitative detection of human ANGPT1 concentrations.

Reactivity Human

Range 31.2 pg/ml - 2000 pg/ml

Tested Applications ELISA

Target ANGPT1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Serum, Plasma


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID [Q15389](#)

Sensitivity 4 pg/ml

Entrez [284](#)

Expiration Date Please enquire.



Representative
Standard
Curve