


Product Datasheet

Human ANGPTL4 ELISA Kit (orb390827)

Description	Sandwich ELISA for quantitative detection of human ANGPTL4 concentrations.
Reactivity	Human
Range	1250 pg/ml - 80000 pg/ml
Tested Applications	ELISA
Target	ANGPTL4
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9BY76
Sensitivity	600 pg/ml
Entrez	51129
Expiration Date	Please enquire.



Representative
Standard
Curve