


Product Datasheet

Human BMP5 ELISA Kit (orb390835)

Description	Sandwich ELISA for quantitative detection of human BMP5 concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	BMP5
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Bone Tissue Homogenates
Sample Types	Bone Tissue Homogenates, Cell Culture Supernates, Serum
Uniprot ID	P22003
Sensitivity	10 pg/ml
Entrez	653
Expiration Date	Please enquire.



Representative
Standard
Curve