


Product Datasheet

Human BMP7 ELISA Kit (orb390836)

Description	Sandwich ELISA for quantitative detection of human BMP7 concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	BMP7
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Bone Tissue Homogenates, Cell Culture Supernatants, Serum, Plasma, Urine
Sample Types	Bone Tissue Homogenates, Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P18075
Sensitivity	30 pg/ml
Entrez	655
Expiration Date	Please enquire.



Representative
Standard
Curve