

Product Datasheet

Human MIP3 beta ELISA Kit (orb390846)

Description Sandwich ELISA for quantitative detection of human MIP3 beta concentrations.

Reactivity Human

Range 15.6 pg/ml - 1000 pg/ml

Tested Applications ELISA

Target CCL19

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Serum, Plasma


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID [Q99731](#)

Sensitivity 8 pg/ml

Entrez [6363](#)

Expiration Date Please enquire.



Representative
Standard
Curve