

Product Datasheet

Human RANTES ELISA Kit (orb390851)

Description Sandwich ELISA for quantitative detection of human RANTES concentrations.

Reactivity Human

Range 7.8 pg/ml - 500 pg/ml

Tested Applications ELISA

Target CCL5

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates


Sample Types Cell Culture Supernates, Serum, Plasma, Tissue Homogenates

Uniprot ID **P13501**

Sensitivity 4 pg/ml

Entrez **6352**

Expiration Date Please enquire.



Representative
Standard
Curve