


Product Datasheet

Human CRISP3 ELISA Kit (orb390876)

Description	Sandwich ELISA for quantitative detection of human CRISP3 concentrations.
Reactivity	Human
Range	1560 pg/ml - 100000 pg/ml
Tested Applications	ELISA
Target	CRISP3
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Urine, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates, Urine
Uniprot ID	P54108
Sensitivity	50 pg/ml
Entrez	10321
Expiration Date	Please enquire.



Representative
Standard
Curve