


Product Datasheet

Human Interferon gamma ELISA Kit (orb390886)

Description	Sandwich ELISA for quantitative detection of human Interferon gamma concentrations.
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	CXCL9
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q07325
Sensitivity	2 pg/ml
Entrez	4283
Expiration Date	Please enquire.



Representative
Standard
Curve