

Product Datasheet

Human COL18A1 ELISA Kit (orb390898)

Description Sandwich ELISA for quantitative detection of human COL18A1 concentrations.

Reactivity Human

Range 78 pg/ml - 5000 pg/ml

Tested Applications ELISA

Target COL18A1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates


Sample Types Cell Culture Supernates, Serum, Plasma, Tissue Homogenates

Uniprot ID **P39060**

Sensitivity 30 pg/ml

Entrez **80781**

Expiration Date Please enquire.



Representative
Standard
Curve