


Product Datasheet

Human AHSG ELISA Kit (orb390899)

Description	Sandwich ELISA for quantitative detection of human AHSG concentrations.
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	AHSG
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P02765
Sensitivity	15 pg/ml
Entrez	197
Expiration Date	Please enquire.



Representative
Standard
Curve