


Product Datasheet

Human FGF19 ELISA Kit (orb390900)

Description	Sandwich ELISA for quantitative detection of human FGF19 concentrations.
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	FGF19
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	O95750
Sensitivity	10 pg/ml
Entrez	9965
Expiration Date	Please enquire.



Representative
Standard
Curve