

Product Datasheet

Human FLT3LG ELISA Kit (orb390908)

Description Sandwich ELISA for quantitative detection of human FLT3LG concentrations.

Reactivity Human

Range 15.6 pg/ml - 1000 pg/ml

Tested Applications ELISA

Target FLT3LG

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Serum, Plasma


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **P49771**

Sensitivity 10 pg/ml

Entrez **2323**

Expiration Date Please enquire.



Representative
Standard
Curve