


## Product Datasheet

### Human IGFBP2 ELISA Kit (orb390925)

|                            |                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Sandwich ELISA for quantitative detection of human IGFBP2 concentrations.                       |
| <b>Reactivity</b>          | Human                                                                                           |
| <b>Range</b>               | 62.5 pg/ml - 4000 pg/ml                                                                         |
| <b>Tested Applications</b> | ELISA                                                                                           |
| <b>Target</b>              | IGFBP2                                                                                          |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.   |
| <b>Note</b>                | For research use only                                                                           |
| <b>Application notes</b>   | Detection Method: Colorimetric<br>Sample Types: Cell Culture Supernatants, Serum, Plasma, Urine |
| <b>Sample Types</b>        | Cell Culture Supernates, Serum, Plasma, Urine                                                   |
| <b>Uniprot ID</b>          | <b>P18065</b>                                                                                   |
| <b>Sensitivity</b>         | 10 pg/ml                                                                                        |
| <b>Entrez</b>              | <b>3485</b>                                                                                     |
| <b>Expiration Date</b>     | Please enquire.                                                                                 |



Representative  
Standard  
Curve