


Product Datasheet

Human IGFBP7 ELISA Kit (orb390928)

Description	Sandwich ELISA for quantitative detection of human IGFBP7 concentrations.
Reactivity	Human
Range	78 pg/ml - 5000 pg/ml
Tested Applications	ELISA
Target	IGFBP7
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	Q16270
Sensitivity	40 pg/ml
Entrez	3490
Expiration Date	Please enquire.



Representative
Standard
Curve