

Product Datasheet

Human IL1 alpha ELISA Kit (orb390930)

Description Sandwich ELISA for quantitative detection of human IL1 alpha concentrations.

Reactivity Human

Range 3.9 pg/ml - 250 pg/ml

Tested Applications ELISA

Target IL1A

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates


Sample Types Cell Culture Supernates, Serum, Plasma, Tissue Homogenates

Uniprot ID **P01583**

Sensitivity 2 pg/ml

Entrez **3552**

Expiration Date Please enquire.



Representative
Standard
Curve