


Product Datasheet

Human IL17RC ELISA Kit (orb390933)

Description	Sandwich ELISA for quantitative detection of human IL17RC concentrations.
Reactivity	Human
Range	250 pg/ml - 16000 pg/ml
Tested Applications	ELISA
Target	IL17RC
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q8NAC3
Sensitivity	15 pg/ml
Entrez	84818
Expiration Date	Please enquire.



Representative
Standard
Curve