


Product Datasheet

Human IL22 ELISA Kit (orb390938)

Description	Sandwich ELISA for quantitative detection of human IL22 concentrations.
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL22
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	Q9GZX6
Sensitivity	15 pg/ml
Entrez	50616
Expiration Date	Please enquire.



Representative
Standard
Curve