


Product Datasheet

Human KLK14 ELISA Kit (orb390945)

Description	Sandwich ELISA for quantitative detection of human KLK14 concentrations.
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	KLK14
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Cell Lysates, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	Q9P0G3
Sensitivity	10 pg/ml
Entrez	43847
Expiration Date	Please enquire.



Representative
Standard
Curve