

Product Datasheet

Human HAVCR1 ELISA Kit (orb390947)

Description Sandwich ELISA for quantitative detection of human HAVCR1 concentrations.

Reactivity Human

Range 15.6 pg/ml - 1000 pg/ml

Tested Applications ELISA

Target HAVCR1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Cell Lysates, Serum, Plasma, Urine


Sample Types Cell Culture Supernates, Serum, Plasma, Cell Lysates, Urine

Uniprot ID [Q96D42](#)

Sensitivity 7 pg/ml

Entrez [26762](#)

Expiration Date Please enquire.



Representative
Standard
Curve