


Product Datasheet

Human LUM ELISA Kit (orb390952)

Description	Sandwich ELISA for quantitative detection of human LUM concentrations.
Reactivity	Human
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	LUM
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P51884
Sensitivity	10 pg/ml
Entrez	4060
Expiration Date	Please enquire.



Representative
Standard
Curve